

# Determination of Alendronate in Human Urine by High-Performance Liquid Chromatography with fluorescence Detection

*Seung Kown Yang<sup>1</sup>, Seungwoo Kang<sup>1,3</sup>, Young Chul Shin<sup>1</sup>,  
Sang Beom Han<sup>3</sup>, Kyung Ryul Lee<sup>1,2</sup>, Hee Joo Lee<sup>1,2</sup>*

*<sup>1</sup>Department of Drug Development Service, BioCore Co., Ltd,*

*<sup>2</sup>Department of Pharmacokinetics, Seoul Medical Science Institute,  
Seoul Clinical Laboratories (SCL)*

*<sup>3</sup>College of Pharmacy, Chung-Ang University*

# Introduction


A fast and sensitive liquid chromatographic method for the determination of alendronate in human urine was developed and validation. Alendronate and internal standard in human urine were extracted by diethylamine solid-phase extraction, 9-fluorenylmethyl derivative. The column effluent was monitored by fluorescence detection at an excitation 260nm and emission wavelengths 310nm. Analytes were separated on a **Onyx C18 (4.6 X 100mm)** with 25mM Sodium pyrophosphate in 20 mM citric acid (pH 3.88) / ACN / MeOH = 72.5 / 23.5 / 4 (v/v/v), as mobile phase. The flow rate was 2 ~ 3 mL/min with the total analysis time of 7.5 min. The standard calibration curves were linear over the concentration range 10-2000 ng/mL with correlation coefficient of 0.999. The limit of quantitation (at signal-to-noise ratio S/N=10) was 10 ng/mL. This method has good precision (intra-day CV(%)  $\leq$  6.47, inter-day CV(%)  $\leq$  5.65) and accuracy (90.56-105.30%) over a wide dynamic range (10-2000 ng/mL). This method has been successfully applied to the pharmacokinetic study of alendronate in human urine.

# Summary of The Analytical Method



---

- Matrix : Human urine
- Sample preparation : SPE, 9-fluorenylmethyl derivative
- Concentration range : 10 to 2000 ng/mL
- Chromatography : HPLC
- Detection mode : FLD
- Quantitation method : Internal standard method
- Quantitation by : Peak area ratio

# 9-Fluorenylmethyl Derivatives



**basic condition**





# Monolithic HPLC Column

## What is Onyx™?

**ONYX**  
Finish First™

### Bimodal Pore Structure


Onyx™ is a silica-based monolithic HPLC column. This technology creates highly porous rods of silica with a revolutionary bimodal pore structure.

The single piece of high-purity polymeric silica gel is then clad in PEEK tubing to make the finished product.

#### Macroporous Structure

**Allows rapid flow (up to 9mL/min) at low pressures**

Each macropore is on average 2  $\mu\text{m}$  in diameter and together form a dense network of pores through which the mobile phase can rapidly flow at low pressure dramatically reducing separation time.



#### Mesoporous Structure

**Creates large surface area**

The mesopores form the fine porous structure (130Å) of the column interior and create a very large surface area on which adsorption of the target compounds can occur.

The unique combination of macropores and mesopores enables Onyx™ monolithic HPLC columns to provide excellent separations in a fraction of the time compared to a standard particulate column.



# A Comparison between Monolithic and ODS(I)


Performance obtained for alendronate using silica-based reversed-phase column and monolithic columns

	$k$	$A_s$	$N_{df}$	$SIN^a$
Conventional silica-based reversed-phase column(ODS)	4.91	1.17	10847.31	12
Onyx monolithic column	3.33	1.03	3424.84	38


# Chromatographic Conditions of HPLC

Mobile phase	25mM Sodium pyrophosphate in 20 mM citric acid (pH 3.88) / ACN / MeOH = 72.5 / 23.5 / 4 (v/v/v)
Flow rate	2 ~ 3 mL/min
Detection	Ex : 260 nm, Em : 310 nm
Injection volume	5 $\mu$ L
Oven temperature	20 $^{\circ}$ C
Column	Onyx C18 (4.6 X 100mm)

# A Comparison between Monolithic and ODS(II)



(a) Chromatogram of human urine using a monolithic column




(b) Chromatogram of human urine using an ODS column




# Calibration Curve of Alendronate

Alendronate




# HPLC Chromatograms of Urine (I)




(a) Chromatogram of blank human urine

# HPLC Chromatograms of Urine (II)



(b) Chromatogram of human urine spike with alendronate (1000ng/mL) and Internal standard


# HPLC Chromatograms of Urine (III)



(c) Chromatogram of human urine sample after oral administration of 70mg alendronate capsule



# Urinary Excretion Rate of Alendronate in Human Urine- Time Curve



# **Conclusion**

---

**We developed a fast and sensitive HPLC-FLD method for the determination of alendronate in human urine. Validation experiments have shown that the assay has good precision and accuracy over a wide concentration range (10-2000 ng/mL). This method is accurate, reproducible and suitable for the analysis of alendronate in clinical samples.**