

# SPE VACUUM MANIFOLDS

- Multiple sample processing saves time and money
- Female Luer inlets fit all male Luer tipped SPE tubes and cartridges
- Inert Teflon® and polypropylene flow path
- Glass chamber allows easy visual monitoring
- Adjustable racks accommodate a variety of test tubes, vials and flasks



## ORDERING INFORMATION

### 10 - Position Tall-Boy Vacuum Manifold\*

Part No.	Description	Unit	Price
AH0-7502	SPE 10-Position Tall-Boy Vacuum Manifold, complete assembly	ea	
AH0-7503	SPE 10-Position Tall-Boy Vacuum Manifold, Glass Chamber	ea	
AH0-7504	SPE 10-Position Tall-Boy Vacuum Manifold, Cover, Gasket and 10 Stopcocks	ea	

### 12 - Position Vacuum Manifold\*

AH0-6023	SPE 12-Position Vacuum Manifold Set, complete assembly	ea	
----------	--	----	--

### 12 - Position Vacuum Manifold Replacement Parts

AH0-6025	SPE 12-Position Glass Chamber	ea	
AH0-6027	SPE Cover, Gasket and 12 Stopcocks	ea	
AH0-6029	SPE Gaskets	2/pk	
AH0-6031	SPE Vacuum Gauge, Valve and Glass Chamber	ea	
AH0-6037	SPE Collection Rack, including shelves, legs, clips and posts	ea	
AH0-6040	SPE Plate for 13 mm Test Tube	ea	
AH0-6043	SPE Plate for 16 mm Test Tube	ea	
AH0-6042	SPE Plate for Volumetric Flask	ea	
AH0-6045	SPE Plate for Autosampler Vial	ea	
AH0-6046	SPE Plate, Dimple	ea	
AH0-6047	SPE Plate, Base	ea	
AH0-6048	SPE Stopcocks	12/pk	
AH0-6052	SPE 12-Position Vacuum Waste Container, polypropylene	10/pk	

### 24 - Position Vacuum Manifold\*

AH0-6024	SPE 24-Position Vacuum Manifold Set, complete assembly	ea	
----------	--	----	--

### 24 - Position Vacuum Manifold Replacement Parts


AH0-6026	SPE Glass Chamber	ea	
AH0-6028	SPE Cover, Gasket and 24 Stopcocks	ea	
AH0-6030	SPE Gaskets	2/pk	
AH0-6032	SPE Vacuum Gauge	ea	
AH0-6038	SPE Collection Rack, including shelves, legs, clips and posts	ea	
AH0-6041	SPE Plate for 13 mm Test Tube	ea	
AH0-6044	SPE Plate for 16 mm Test Tube	ea	
AH0-6049	SPE Stopcocks	24/pk	

### General Vacuum Manifold Accessories

AH0-6033	SPE Manifold Needles, polypropylene	12/pk	
AH0-6034	SPE Manifold Needles, polypropylene	24/pk	
AH0-6035	SPE Manifold Needles, stainless steel	12/pk	
AH0-6036	SPE Manifold Needles, stainless steel	24/pk	
AH0-6050	SPE Drying Attachment for 12-Position manifold	ea	
AH0-6051	SPE Drying Attachment for 24-Position manifold	ea	
AH0-6053	Female Luer Fittings	2/pk	
AH0-6054	Male Luer Fittings	2/pk	
AH0-6055	Support Posts for Rack	3/pk	
AH0-6056	Legs for Cover, black	4/pk	
AH0-6057	Vacuum Gauge and Valve Assembly	ea	
AH0-6058	Valve Assembly Only	ea	
AH0-6059	Vacuum Gauge Only	ea	
AH0-6060	Retaining Clips	12/pk	
AH0-6061	Vacuum Manifold Plugs	50/pk	
AH0-6062	Control Valve, Teflon®	25/pk	
AH0-6063	Control Valve, Teflon®	50/pk	
AH0-6064	Teflon® Needles	100/pk	
AH0-6065	Teflon® Needles	500/pk	
AH0-7191	Adaptor Caps for 1, 3 and 6 mL SPE tubes, polyethylene, with Luer hole in the top	15/pk	
AH0-7378	Adaptor Caps for 12 and 20 mL SPE tubes, polyethylene, with Luer hole in the top	5/pk	
AH0-7379	Adaptor Caps for 60 mL SPE tubes, polyethylene, with Luer hole in the top	5/pk	
AH0-7907	Adaptor Caps for 150 mL SPE tubes, polyethylene, with Luer hole in the top	ea	

**Have Large Sample Volume but Need a Small Bed Mass?**

**Use an adaptor cap to attach another SPE tube, which can be used to increase the reservoir size for washing or eluting solvents.**



See p. 16 for Strata Impact Precipitation Plates.

\* Manifolds include: Vacuum-tight glass chamber, polypropylene lid with gasket, bleed valve and gauge, stopcock valves, collections racks, polypropylene needles



(1) The 10-position Tall Boy Vacuum Manifold Collection Rack includes 4 plates: one base plate, one dimple plate, one small plate and one large plate and three riser bar legs, along with 12 manifold clips to support the plates. The assembly also includes 10 polypropylene needles, 10 stopcocks and 4 black legs to support the lid when taken off the glass block.

(2) The 12-position Collection Rack consists of 3 support posts, bottom plate, 13 mm plate, 16 mm plate, autosampler plate, volumetric plate, and 12 retaining clips. (3) The 24-position Collection Rack consists of 3 support posts, bottom plate, dimple plate, 13 mm plate, 16 mm plate, and 24 retaining clips.